



**Electric Delivery Bed
HW502A**

- Used for Obstetrics and Gynecology Delivery, Gynecological Examination, Diagnosis, Surgery, Emergency Caesarean Section
- Hand controller or foot controller
- Built-in auxiliary surgery plate

Technical Specification

Table size	1850mm*650mm
Table height movement	600-900mm
Trendelenburg/Reverse trendelenburg	20°/35°
Back plate adjustment	0°-85°
Auxiliary plate size	550mm*580mm



**Electric Delivery Bed
HW502A1**

- Equip with foot controller
- Optional examination light which can be used for gynecology examination.

Technical Specification

Table size	1900mm*600mm
Table height movement	600-900mm
Back plate adjustment	70°/-10°
Trendelenburg/Reverse trendelenburg	20°/20°
Foot plate horizontal movement	90°
Load capacity	250kgs



LDR Bed
HW502C

- LNAK linear motor
- Central braking system
- Guardrail operation buttons
- CPR Cardiac Resuscitation Rescue Function(back plate drop rapidly)

Technical Specification

Table size	2200mm*1100mm
Mattress size	980-900mm
Table height movement	480mm-845mm
Guardrail height movement	405mm
Trendelenburg	8°
Back plate adjustment	65°
Seat plate adjustment	12°
Leg plate height movement	170mm
Leg plate outward and up	90°



LDR Bed
HW502C1

- Reliable performance and stable low noise
- Hand controller and foot controller
- Sound system for option

Technical Specification

Table size	2050mm*980mm
Table height movement	580-950mm
Back plate adjustment	70°
Trendelenburg/Reverse trendelenburg	5°/12°
Foot plate height movement	0-550mm
Diameter of wheel	125mm