

Electro-hydraulic Operating Table HW503E

Technical Specification

Table top size	2100mm*520mm
Table height	750-1100mm
Trendelenburg/Reverse Trendelenburg	25°/25°
Lateral Tilt(L/R)	20°/20°
Back Plate(up/down)	80°/25°
Longitudinal slide	400mm
Head Plate(up/down)	90°/-90°
Leg Plate(up/down)	15°/-90°
Built-in kidney bridge	120mm





Electric Operating Table HW503D

- German 12cm wheels
- Auxilary control board
- Memory foam mattress
- Back up battery

Technical Specification	
Table top size	2100mm*550mm
Table height	710-1010mm
Longitudinal slide	350mm
Trendelenburg/Reverse Trendelenburg	25°/25°
Lateral Tilt(L/R)	15°/15°
Head Plate(up/down)	40°/-90°
Leg Plate(up/down)	15°/-90°
Back Plate(up/down)	80°/-20°
Built-in kidney bridge	110mm



Electric Operating Table HW503B

- Two configurations for option
- One configuration has longitudinal slide function which can be used with C-arm
- One configuration has lower height which can be used with neurosurgery

Technical Specification		
Mode	HW-503B	HW-503B1
Table top size	2000mm*500mm	2000mm*500mm
Height	750-1000mm	550-800mm
Built-in kidney bridge	100mm	80mm
Trendelenburg/Reverse Trendelenburg	25°/25°	25°/20°
Lateral Tilt(L/R)	20°/20°	20°/20°
Back Plate(up)	70°/-35°	45°/-20°
Leg Plate(up/down)	15°/-90°	15°/-90°
Head Plate(up/down)	45°/-90°	30°/-90°
Longitudinal slide	350mm	



Electro-hydraulic Operating Table HW503A

Low height especially for neurosurgery

- Built-in kidney bridge
- Optional longitudinal slide

Technical Specification

Table top size	2000mm*520mm
Table height	600-900mm
Longitudinal slide	400mm
Trendelenburg/Reverse Trendelenburg	22°/22°
Lateral Tilt(L/R)	20°/20°
Head Plate(up/down)	90°/-90°
Leg Plate(up/down)	15°/-90°
Back Plate(up/down)	80°/-20°



Image integrated Operating Table HW503E1

- 35mm thickness carbon fiber tabletop
- Advanced electric-mechanic transmission technology
- Longitudinal slide 600mm to 800mm

Technical Specification

Table top length	2600mm
Table top width	550mm
Longitudinal slide	600mm-800mm
Horizontal slide	110mm